



张家港麦迪卡医疗器械有限公司  
ZHANGJIAGANG MEDI MEDICAL EQUIPMENT CO., LTD.



CE ISO9001 ISO13485



**Model: YA-XD1C Electric Ophthalmic operation table**

**Dimension:** 1950mm\*520mm\*550-800mm

**Functions:**

- 1.HI-LOW: 550mm-800mm
- 2.Movement of head section  $\geq 100$ mm

**Structure:** Stainless steel

**Accessories:**

- 1 pair armrests for patient, 1 pair armrests for doctor  
1pc anesthesia screen, 1 pair belts for fixing armrests

Maximum loading capacity: 180kgs

Packing:1.1M3                      G.W.: 238KGS